# wellion®

## WAVE *plus*

# HANDBUCH

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# OWNERS GUIDE

<b>TABLE OF CONTENTS</b>	<b>25</b>
<b>ABOUT BLOOD PRESSURE</b>	<b>26</b>
What is blood pressure?	26
What is hypertension and how is it controlled?	26
Why measure blood pressure at home?	26
WHO blood pressure classification	26
Blood pressure variations	27
<b>NOTES ON SAFETY</b>	<b>27</b>
Explanation of the icons	27
<b>PRECAUTIONS BEFORE USE</b>	<b>29</b>
<b>PARTS IDENTIFICATION</b>	<b>29</b>
<b>INSERT OR REPLACE BATTERIES</b>	<b>30</b>
Symbols on display	30
Accessory	30
<b>IRREGULAR HEARTBEAT SIGNAL (♥)</b>	<b>31</b>
<b>TIME AND SYSTEM SETUP</b>	<b>32</b>
<b>SETUP UNIT mmHg/kPa</b>	<b>32</b>
<b>WHO BLOOD PRESSURE CLASSIFICATION DISPLAY</b>	<b>33</b>
<b>ATTACHING THE ARM CUFF</b>	<b>33</b>
How to take proper measurements	34
<b>HOW TO MEASURE BLOOD PRESSURE</b>	<b>35</b>
Read memory	35
Delete Memory	35
<b>CARE AND MAINTENANCE</b>	<b>36</b>
<b>TROUBLESHOOTING</b>	<b>37</b>

SPECIFICATIONS	38
STATEMENT	38
FEATURE OF THE PRODUCT	39
REQUESTS FROM MANUFACTURER	39
REGARDING EMC STATEMENT	40
MANUFACTURER´S WARRANTY	40

## **INTRODUCTION**

The Wellion WAVE PLUS Blood Pressure Monitor uses the oscillometric method of blood pressure measurement.

The Measurement Automatic Electronic Blood Pressure Monitor is intended for use by medical professionals or at home to monitor and display diastolic, systolic blood pressure and pulse rate, with an arm cuff around the left upper arm according to the instruction in the "ATTACHING THE ARM CUFF". The expected life of the product is 5 years.

The product complies with the electromagnetic compatibility requirements of EN 60601-1-2 and safety standards of EN 60601-1 and performance of EN 80601-2-30 as specified in EEC directive 93/42/EEC.

## ABOUT BLOOD PRESSURE

### 1. What is blood pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimetres of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

### 2. What is hypertension and how is it controlled?

Hypertension is an abnormally high arterial blood pressure. If this is left unattended it can cause many health problems including stroke and heart attack.

Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.

To prevent hypertension or to keep it under control:


- Do not smoke
- Exercise regularly
- Reduce salt and fat intake
- Have regular physical checkups
- Maintain proper weight

### 3. Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than the value measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

### 4. WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), and shown in chart beside.

## 5. Blood pressure variations

An individual's blood pressure varies on daily and seasonal basis. It may vary by 30 to 50 mmHg due to various conditions during the day.


In hypertensive individuals, variations are even more pronounced.

Normally, the blood pressure rises during work or while playing and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement.

Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.

Be sure to note date and time when recording your blood pressure.

Consult your doctor to interpret your blood pressure data.

### NOTES ON SAFETY

The warning signs and sample icons shown here are listed for your safety and the correct use of the unit, so as to prevent injuries or damages to the device.

The icons and meanings are as follows:

#### Explanation of the icons

 The icon indicates prohibitions (what you should not do). Matters involving actual prohibitions are indicated by text or pictures in or near. The left icon refers to “general prohibition”.

 The icon indicates something that is compulsory (what must always be observed). Matters involving actual mandatory actions are indicated by text or pictures in or near “”. The left icon refers to “general compulsion”.

 The icon indicates something can't be disassembled or “Don't disassemble”.

 Type BF Applied part

Attention and read before use

Please refer to the instructions for use

Product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling. Comply with WEEE Directive 2012/19/EU.

Indicates a medical device that needs to be protected from moisture.

**MR Unsafe**  
Do not use this  
equipment in the MRI  
scan room



<p>Patient must follow doctor's instruction and should not perform self-judgment and self-treatment by the measuring result. Self-diagnosis of measured results and treatment are dangerous. The device should not be used to judge illness, first aid measures and continuous monitoring of blood pressure. This device cannot be used for measurements on patients who are under transport or under surgical treatment.</p> <p>Please press the "ON/OFF" button to stop the device from working when you feel uncomfortable with the arm, or if the air is inflating abnormally without stop.</p>	
<p>Do not let a child below the age of 12 years and people who can't express one's intention use the device. When it is used by people at the age of 12-18 years, an adult should accompany it, as it may cause an accident or trouble.</p>	
<p>Do not use the device for purposes other than measuring blood pressure, as it may cause an accident or trouble.</p>	
<p>Please do not use mobile phones around the device. Please do not use the device around a magnetic field.</p>	
<p>The device is prohibited from being used during movement.</p>	
<p>Do not use the equipment in outdoor or shower rooms.</p>	
<p>Stop using the equipment immediately, if it is in contact with water.</p>	
<p>Do not disassemble, repair, or remodel the main unit or the cuff of the blood pressure monitor. This will cause the unit to function erroneously.</p>	

## PRECAUTIONS BEFORE USE


1. If you are taking medication, consult with your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your doctor.
2. For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be fluctuation in blood pressure values measured at the upper arm versus at the wrist.
3. Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields.  
To prevent such interference, use the monitor at a sufficient distance from such devices or turn them off.
4. Before using you should wash your hands.
5. Do not measure on the same arm which uses simultaneously another medical monitoring ME Equipment, otherwise it could cause loss of function.
6. Consult your doctor if unexpected readings are obtained, also please refer to "TROUBLESHOOTING" of the manual.
7. The reading at home is probably a little lower than measured in the hospital due to the steady mood at home.
8. Cuff pressure range 0-280 mmHg

## PARTS IDENTIFICATION


LCD Display  
(Date and Time,  
Systolic Blood Pressure,  
Diastolic Blood Pressure,  
Pulse/min)


### Symbols on display

### Accessory

Manual


AC Adapter (optional)

Use exclusively a safety-proofed AC adapter which can be purchased from the manufacturer of the blood pressure meter.

### INSERT OR REPLACE BATTERIES

1. Remove the battery cover.
2. Insert new batteries into the battery compartment as shown, taking care that the polarities (+) and (-) are correct.
3. Close the battery cover. Use only LR6, AA batteries.

Disposal of empty battery to the authorized collecting party subject to the regulation of each individual territory.

 Caution:

- Insert the batteries as shown in the battery compartment. If not, the device will not work.
- When „LOW BATTERY“ symbol  blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
- „LOW BATTERY“ symbol  does not appear when the batteries run out.
- Please ensure to distinguish positive polar (+) and negative polar (-) of batteries when replacing batteries.
- Batteries, which have a fluid on the surface or have been modified, must not be inserted into the products.
- Battery short circuit must be prevented.
- Battery life varies with the ambient temperature and may be shorten at low temperatures.
- The batteries may leak and cause a malfunction.
- Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a shorter life.
- Used batteries may leak and damage the main unit. Please observe the following points:
  - If you are not going to use the unit for a long period of time (approximately three months or more), remove the batteries.
  - Replace empty batteries with their polarities in the correct direction.


**IRREGULAR HEARTBEAT SIGNAL «♥»**

When the Wellion WAVE PLUS detects an irregular heartbeat rhythm, two or more times during the measurement, the irregular heartbeat symbol «♥» will appear on the display with the measurement value.

An irregular heartbeat rhythm is defined as a rhythm that is 25% less or 25% more than the average rhythm detected while the monitor is measuring the systolic and diastolic blood pressure.

## TIME AND SYSTEM SETUP


1. Press "SET" button for the Time display.
2. In the off state, press and hold the "SET" button until the year number displays and flashes on the LCD to enter setting mode.
3. Press "MEM" button to adjust the year, then press "SET" button again to save your setting and enter the month-setting mode.
4. Press "MEM" button to adjust the month. Repeat procedure to adjust date/hour/minute.

## SETUP UNIT mmHg/kPa

The results can be displayed as one of two different units: mm Hg (mmHg) or kPa (kPa).

The preset factory unit is mmHg. Press "ON / OFF" button for 10 seconds to display unit switching interface, then press "MEM" button to select mmHg / kPa, press "ON / OFF" button to exit. After the units have been selected and confirmed according to the instructions above, the display shows mmHg / kPa. After a normal start, the selected units are then shown as blood pressure units. Please note that all results`units in the device memory will be changed accordingly as well.

## WHO BLOOD PRESSURE CLASSIFICATION DISPLAY


Diastolic blood pressure Reference material:  
Journal of hypertension 1999. vol 17 No.2

■	Grade 3 hypertension (severe)
■	Grade 2 hypertension (moderate)
■	Grade 1 hypertension (mild)
■	High-normal
■	Normal
■	Optimal

EN

### ATTACHING THE ARM CUFF


1. Wrap the arm cuff around the upper arm, about 2-3 cm above the elbow, as shown. Place the cuff direct on the skin, as clothing may cause a faint pulse and result in a measurement error.
2. Constriction of the upper arm, caused by rolling up a shirtsleeve, may prevent accurate readings.
3. Secure the arm cuff with the Velcro Strip in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the arm cuff is at the same height as the heart. Make sure that the tube is not kinked.
4. Measure your arm circumference for the cuff selection, refer to "SPECIFICATIONS".

## How to take proper measurements


For best accuracy in blood pressure measurement:

- Patient in seated position, including
  - 1) comfortably seated
  - 2) legs uncrossed
  - 3) feet flat on the floor.
  - 4) back and arm supported
  - 5) middle of the CUFF at the same level as the heart
- Remain still and keep quiet during measurement.
- Relax as much as possible and do not talk during the measurement process.
- Measure your blood pressure at about the same time every day.
- Do not measure right after physical exercise or after a bath. Take a rest for twenty or thirty minutes before taking the measurement.
- Other conditions which could affect the measurement results:
  - Within in an hour after dinner, after having wine ,coffee, red tea, sports, bathing; talking, being nervous, being in unsteady mood, bending forward, moving, room temperature dramatically changing during measuring; In the moving vehicles, long time continuous measuring.


## HOW TO MEASURE BLOOD PRESSURE

1. Set up the arm cuff to your upper arm as described in the previous section of "ATTACHING THE ARM CUFF".
2. Press the "ON/OFF" button, all icons automatically appear for two seconds on the display, then the device switches to measurement and displays "00".

3. When the measurement starts, the cuff in the strap will automatically inflate. This is indicated by the "Inflation" symbol "♥". When the measurement is completed, the results will be displayed on the LCD display.

### Read memory

Press "MEM" button to inquire memory average values "AU9" Average Value Display: the latest 3 groups of memory average values (Memory values are displayed regardless of period). Press the "MEM" button to scroll through the memorized results, starting with the latest measurement. Press "MEM" to go up (latest result) and "SET" to go down (oldest) in the list of results.

### Delete Memory


Enter the "READ MEMORY" mode, press and hold the "MEM" button for five seconds until the LCD displays "no". This indicates that all stored results have been deleted successful.

## CARE AND MAINTENANCE

To keep your digital blood pressure monitor in the best condition and to protect the unit from damage, follow the directions listed below:

- Keep the monitor in the storage case when not in use.
- Do not fold the arm cuff too tightly.
- The fabric fastener could touch the inner surface of the arm cuff and damage it.
- Clean the monitor and cuff with a soft dry cloth.
- Do not use any abrasive or volatile cleaners.

### ⚠ Caution:

- Do not submerge the device or any of the components in water.
- Do not subject the monitor to extreme hot or cold temperatures, humidity or direct sunlight.
- Store the device and the components in a clean, safe location.
- Do not subject the monitor to strong shocks, such as dropping the unit on the floor.
- Remove the batteries if the unit will not be used for three months or longer.
- Always replace all the batteries with new ones at the same time.

### Information:

- This product is designed for use over an extended period of time. It is generally recommended that it be inspected and calibrated every two years to ensure proper function and performance.

(\*Pressure calibration is done by EU representatives)

## TROUBLESHOOTING

If you have trouble in using the product please check the following points first:

EN

Error display	Possible cause	How to correct
Nothing is displayed when you push the POWER button or battery icon  flashes	Batteries installed	Insert batteries
	Batteries are empty	Replace batteries
	The polarities of batteries placed wrongly	Insert batteries in the correct polarities
E0: Device can't detect the signal pulse	There may be many clothes between the cuff and the upper arm. Please wear a thin shirt when measuring.	Constriction of the upper arm caused by rolling up a shirtleeve, may prevent accurate readings.
	The cuff may be worn in wrong direction and the tube interface doesn't correspond to the arterial position.	
	If the E0 error code is not caused by the above mentioned reasons the device is defect.	
E1: Can't normally Increase pressure	Check that your arm cuff has no air leakage	Replace arm cuff with a new one
E3: inflate pressure too high		Re-measurement or have the device re-calibrated
E2E4: shaking while measurement	Hand or body shaking while measurement	Keep static and correct gesture to measure again
Battery icon  displayed	Batteries have low power	Replace batteries and repeat the measurement
The systolic or diastolic pressure value is too high	The arm cuff was held lower than your heart level	Keep static and keep correct gesture, and repeat measurement
	The arm cuff was not attached properly	
	You moved your body or spoke during measurement	
The systolic or diastolic pressure value is too low	The arm cuff was held lower than your heart level	Keep static and keep correct gesture, and repeat measurement
	You moved your body or spoke during measurement	

## SPECIFICATIONS

**Measuring Method:** Oscillometric Measurement

**Measuring Range:** Pressure: 30 - 280 mmHg  
Pulse: 40 - 199 beats/minute

**Memory:** 90 Memories

**Operating condition:** + 5°C - + 40°C;  
15% RH - 93% RH;  
Atmospheric pressure: 70kPa - 106kPa

**Dimensions:** Approx: 140(W) X 110(H) X 48(D) mm

**Classification:** Type BF

**AC adapter:** Input: 100 – 240V 50/60Hz, 0,2A  
Output: 6V, 500 mA

**Indication:**

**Accuracy:**

**Power supply:**

**Storage condition:**

**Weight:**

**Upper arm circumference:** 22 - 42 cm

Use exclusively a safety-proofed AC adapter which can be purchased from the manufacturer of the blood pressure meter.

Digital LCD display

Static Pressure:  $\pm 3$ mmHg

Pulse:  $\pm 5\%$

4x 1,5V batteries LR6 or AA

**Use alkaline batteries:**

allow over 200 measurements

- 20°C - + 55°C;

0% RH - 93% RH;

Atmospheric pressure: 50kPa - 106kPa

Approx: 410 g, excl. batteries

Note that the specifications may be changed without notice in case of implemented product improvements.

1. Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
2. Degree of protection against electric shock: TYPE BF APPLIED PART.
3. Mode of operation: CONTINUOUS OPERATION.
4. Equipment not suitable for category AP & APG equipment use in presence.

## STATEMENT

The system might not meet its performance specifications if stored or used outside the temperature and humidity as mentioned below:

- Operating conditions: +5°C - +40°C                      15%RH - 93%RH;                      70kPa - 106kPa
- Storage conditions: -20°C - +55°C                      0%RH - 93%RH

## FEATURE OF THE PRODUCT

1. Memory can store 90 measurements.
2. Large and clear LCD display.
3. WHO blood pressure classification display.
4. Easy to use: Press a button to automatically measure; Record the measurement values and measurement time.
5. Automatically turns off (within 1 minute) to save power.

## REQUESTS FROM MANUFACTURER

- Make sure that the connection tube is not kinked before starting the measurement to avoid any injury to patient.
- For any patient, do not measure more than 3 times continuously, it should be at least above 5 minutes of interval rest between any two measurements, otherwise will cause exsanguination in the blood vessel.
- Do not measure your blood pressure over 6 times each day.
- Do not apply the cuff over a wound as this can cause further injury.
- Do not measure on the arm which is on the side of a mastectomy, otherwise it could cause injury.
- Observe the air pressure value from the LCD display.
- When measuring, the device cannot exceed 280 mmHg, otherwise please press the "ON/OFF" button to stop.
- Do not use force to bend the arm cuff or the air tube.
- Do not knock or drop the main unit.
- Always use the specified accessories in the manual, the use of other parts not approved by the manufacturer may cause faults or injuries.
- For service information, parts list etc., please contact the dealer.
- The patient is an intended operator.
- Not servicing and maintenance while the ME EQUIPMENT is in use.
- The user can maintain the product, the maintenance method is described in the maintenance instructions of manual.
- Stop using the equipment immediately, if it is in contact with water.

## **REGARDING EMC STATEMENT**

1. A statement to the effect that this equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS; and
2. a statement to the effect that wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept distanced away from the equipment. The distance is calculated by the MANUFACTURER from the 800 MHz to 2,5 GHz column of Table 5 or Table 6 of EN 60601-1-2:2015, as appropriate.

## **MANUFACTURER'S WARRANTY**

The producer warrants Wellion WAVE PLUS blood pressure meter against defects in materials and workmanship for a period of two years from the date of purchase. The warranty is lost if the instrument is misused, poorly maintained or is tampered with. Liability under this warranty is limited to the repair of defective parts or - at the discretion of the producer - to the replacement of the instrument. The right to rescind the purchase agreement exists only if the replacement is also found to be defective. Claims other than these will not be dealt with. The warranty is not valid if the damage results from misuse, maltreatment, tampering, human error and the use of extreme force. This warranty is valid only if the date, the stamp and the signature of the dealer are recorded on the warranty card on the date of purchase. The warranty period is not to be extended by any claim made under this warranty.

## Appendix 1 Guidance and Manufacturer Declaration Tables

<b>Guidance and manufacturer's declaration- electromagnetic emission</b>		
The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Model PG-800B11 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Model PG-800B11 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complied	

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment-guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/bursts IEC 61000-4-2	± 2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	± 2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV differential mode line-line	$\pm 0.5$ kV, $\pm 1$ kV differential mode line-line	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30% dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) for 250/300 cycle at 0°	0% UT (100% dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30% dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) for 250/300 cycle at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Model PG-800B11 product name requires continued operation during power mains interruptions, it is recommended that the Model PG-800B11 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a. c. mains voltage prior to application of the test level.			

**Guidance and manufacturer's declaration-electromagnetic immunity**

The Model PG-800B11 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800B11 should assure that it is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF  IEC 6100-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands</p>	<p>6 V</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Models PG-800B11, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$ $d = \left[ \frac{12}{E_1} \right] \sqrt{P} \quad \text{80MHz to 800MHz}$ $d = \left[ \frac{23}{E_1} \right] \sqrt{P} \quad \text{800MHz to 2.7GHz}$

Radiated RF IEC 6100-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b Interference may occur in the vicinity of equipment marked with the following symbol: ((⋯)) ▲
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.			
b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.			
c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-800B11 is used exceeds the applicable RF compliance level above, the Model PG-800B11 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-800B11.			
d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			

### Recommended separation distances between portable and mobile RF communications equipment and the Model PG-800B11

The Model PG-800B11 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model PG-800B11 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model PG-800B11 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz ISM bands $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	150 kHz to 80 MHz ISM bands $d = \left[ \frac{12}{V_2} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{12}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[ \frac{23}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.27
100	12	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

NOTE 3: An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Authorised representative: Lotus NL B.V.
Manufacturer: Shenzhen Pango Medical Electronic Co., Ltd.